




INTRODUCING THE NEW GENERATION PTS

INDEPENDENT INTELLIGENCE



PTS PERSONALIZED
TOURNIQUET
SYSTEMS

SAFETY, PERFORMANCE AND PERSONALIZATION IN A COMPACT SOLUTION

2019

TWO-IN-ONE

Independent instruments provide flexibility and cost savings to facilities with multiple ORs, by allowing twin systems to be separated and used as two single units. For IVRA procedures, simply bring units back together

INTELLIGENT TOUCH

Intuitive user interface simplifies key tourniquet parameters and user training

PERSONALIZED PRESSURE LIMITS

Limits normal maximum pressure to 475 mmHg, and allows selection of extended pressure up to 600 mmHg if required for individual patients

ADVANCED LEAK DETECTION DURING SURGERY

Monitors the tourniquet cuff, tubing and connectors, and safely alerts staff to hazards

SPEED

Reduces set-up and cuff inflation times

INTEGRATED TOURNIQUET CUFF TESTING

Enables cuff testing before or after surgery as recommended in current AORN Standards and Recommended Practices

ADVANCED SAFETY FEATURES

Intuitive time and pressure safety alarms; unique self-test and calibration tests; visible displays of pressure and time

GREATEST BATTERY LIFE

Up to 10 hours of operation with tested cuffs

CONNECTIVITY

Allows integration with OR information networks



Delfi's state-of-the-art PTSii electronic tourniquet systems have advanced personalization safety features and save valuable space in modern surgical settings.

Available in both single and twin system configurations, fully independent twin channels are color-coded for optimal safety during simultaneous bilateral limb surgeries and IVRA (Bier Block) procedures.

The PTS ii Personalized Tourniquet System is used with Delfi's unique line of personalized tourniquet cuffs for pediatric to bariatric patients and is compatible with most standard single port tourniquet cuffs.

PTSii Personalized Tourniquet Systems

Single Channel . Single Port System	9-2200-001
Dual Channel . Single Port System	9-2200-500
Single Channel . Single Port System & Easi-Fit Cuffs Kit	9-2200-450
Dual Channel . Single Port System & Easi-Fit Cuffs Kit	9-2200-400
Roll Stand with Basket	9-2200-550

PTSii Instrument Dimensions

Height	Width	Depth	Weight
180mm	120mm	110mm	1.26kgs
7.0in	4.7in	4.3in	2.75lbs

About Us

Delfi Medical Innovations is a world leader in advancing tourniquet safety and technology. Members of Delfi's team have been exclusively involved in tourniquet development and evaluation for more than 40 years.

Delfi's patented tourniquet instruments and tourniquet cuffs have been designed to help surgeons and operating-room nurses provide personalized safe and effective care for their patients.

Contact Us

For sales enquiries and ordering information contact Customer Service at:

1.800.933.3022 (US & Canada)
+ 604.742.0600 (Elsewhere)
sales@delfimedical.com

Delfi Medical Innovations Inc.
106 - 1099 West 8th Avenue
Vancouver, BC Canada V6H 1C3

To learn more on how our tourniquet systems help provide personalized safe and effective care for your patients visit:

delfimedical.com