

PTS TRANSPORT BAG

COMFORT AND EASY ACCESS IN A CARRY-ON

Delfi's **PTS Transport Bag** is custom designed to hold a single **PTS System** including accessories and cuffs. A separate padded instrument case and cuff pouch provide added protection, while hoses and personal items are conveniently stored in a zippered compartment in the lid.

The low profile handle, removable carrying strap and pass-through luggage pocket make the **PTS Transport Bag** the perfect carry-on for clinicians on-the-go.

FEATURES

- Rugged fabric construction
- Rubber feet
- Removable adjustable carrying strap
- Low profile handle
- Pass-through luggage pocket
- Business card holder
- Separate padded PTS case and cuff pouch

PART NOS.

- 9-2200-300 Transport Bag (with Padded Instrument Case & Cuff Pouch)
- 9-2200-301 Instrument Case, Padded
- 9-2200-302 Cuff Pouch, Padded

DIMENSIONS

- Exterior (L×W×D) 16" x 8.5" x 9.5" in
- Interior (L×W×D) 15.3" x 7.5" x 8" in
- Weight 3 Lbs (empty)
- Weight 10 Lbs (with PTS System & Cuffs)



Contact Us

For sales enquiries and ordering information contact Customer Service at:

1.800.933.3022 (US & Canada)
+ 604.742.0600 (Elsewhere)
sales@delfimedical.com

Delfi Medical Innovations Inc.
106 - 1099 West 8th Avenue
Vancouver, BC Canada V6H 1C3

