

PTS ROLL STAND

PERFORMANCE IN COMPACT DESIGN

The [PTS Roll Stand](#) is a strong and sturdy mobile stand for use in a wide range of healthcare and clinical environments, including operating rooms and rehabilitation clinics. Designed to securely hold up to two separate PTS instruments, plus accessories.

The combination of a low-profile base, all-locking castors, and a solid aluminum column improves stability, and increases maneuverability. The included [PTS Roll Stand Basket](#) stores the tourniquet cuffs, keeps hoses tidy and makes access easy.

FEATURES

- Sturdy, mobile roll stand
- Five locking castors
- Accessory wire basket included, featuring dual hose management hooks
- Small foot print
- Easy to clean design

DIMENSIONS

- 50.25" x 22" (Assembled)
- Weight 14 Lbs

PART NO.

- 9-2200-550



Contact Us

For sales enquiries and ordering information contact Customer Service at:

1.800.933.3022 (US & Canada)
+ 604.742.0600 (Elsewhere)
sales@delfimedical.com

Delfi Medical Innovations Inc.
106 - 1099 West 8th Avenue
Vancouver, BC Canada V6H 1C3

