

MLPS

Matching Limb Protection Sleeves

Delfi's Matching Limb Protection Sleeves (MLPS) are uniquely matched to the size and shape of each specific Delfi cuff model, to provide optimal personalized limb protection.

Delfi's MLPS help protect underlying skin and soft tissues from injuries such as pinching, wrinkling and shearing. The use of Delfi MLPS is supported by published surgical literature.


Every Delfi sterile single-use and reusable cuff includes one Matching Limb Protection Sleeve (MLPS). For MLPS selection refer to Purchasing Guide.



FEATURES

- Quick and easy application
- Soft and wrinkle-free
- Two layers of protection
- Prevents skin entrapment at cuff overlap
- Latex-free
- Save time with colour coding for easy selection
- Protects skin from pinching and shearing
- Does not clog hook fasteners
- Helps extend the life of reusable cuffs
- Non slippery material prevents cuff from sliding

LIMB PROTECTION COMPARISON



As recommended in current AORN Standards, Recommended Practices and Guidelines, using appropriate limb protection under the cuff can significantly reduce soft tissue damage

In extensive testing using a molding and digitizing technique, results have shown that a fitted sleeve that is matched to a specific cuff and that applies a controlled amount of compression to the limb virtually eliminates pinching and wrinkling of the underlying skin reducing localized pressure gradients, which can result in soft tissue damage.

Matching Limb Protection Sleeves

EASI-FIT MLPS

For use with Easi-Fit Tourniquet Cuffs


Part No.	Limb Circumference Range (10 per box)
9-7920-018	18" Cuffs (11" - 16")
9-7920-024	24" Cuffs (15" - 22")
9-7920-034	34" Cuffs (21" - 30")
9-7920-044	44" Cuffs (29" - 40") (5 per box)
9-7920-052	52" Cuffs (38" - 46") (5 per box)



VARI-FIT MLPS

For use with Vari-Fit Tourniquet Cuffs

Part No.	Limb Circumference Range (10 per box)
9-7920-002	Arm (11" - 16")
9-7920-003	Thigh (15" - 22")
9-7920-001	Lower Leg (21" - 30")
9-7920-002	IVRA (arm) (29" - 40")



PETITE-FIT MLPS

For use with Petite-Fit Tourniquet Cuffs

Part No.	Limb Circumference Range (10 per box)
9-7920-002	3" Cuffs (7" - 12 1/2")
9-7920-003	3.5" Cuffs (8 1/4" - 15 1/4")



PEDI-FIT MLPS

For use with Pedi-Fit Tourniquet Cuffs

Part No.	Limb Circumference Range (10 per box)
9-7900-001	1.5" Cuffs (3 1/2" - 5 3/4")
9-7900-002	2.25" Cuffs (5" - 8 1/4")
9-7900-003	3" Cuffs (7" - 12 1/2")
9-7900-004	3.5" Cuffs (8 1/4" - 15 1/4")



About Us

Delfi Medical Innovations is a world leader in advancing tourniquet safety and technology. Members of Delfi's team have been exclusively involved in tourniquet development and evaluation for more than 40 years.

Delfi's patented tourniquet instruments and tourniquet cuffs have been designed to help surgeons and operating-room nurses provide personalized safe and effective care for their patients.

Contact Us

For sales enquiries and ordering information contact Customer Service at:

1.800.933.3022 (US & Canada)
+ 604.742.0600 (Elsewhere)
sales@delfimedical.com

Delfi Medical Innovations Inc.
106 - 1099 West 8th Avenue
Vancouver, BC Canada V6H 1C3

To learn more on how our tourniquet systems help provide personalized safe and effective care for your patients visit:

delfimedical.com